



KEEPING LIFE-SAVING TECHNOLOGY AT THE
Heart of your Business



CINTAS
READY FOR THE WORKDAY[®]



GOOD SAMARITAN LAWS

All 50 states and the District of Columbia have Good Samaritan-related laws. The details of Good Samaritan laws/acts vary by jurisdiction, including who is protected from liability and under what circumstances. Although these laws differ by state, all are designed to encourage people to help a stranger who needs assistance by reducing or eliminating the fear that, if they do so, they will suffer possible legal repercussions if they inadvertently make a mistake in treating the victim.

OSHA RECOMMENDS AEDS IN THE WORKPLACE

"AEDS SAVE LIVES! These devices have a proven track record of saving lives in public places as well as in the workplace. They can do the same for you and your employees. Please consider installing AEDs in your workplace."

- U. S. Department of Labor, Occupational Safety and Health Administration (OSHA), 3185-09N 2003
<https://www.osha.gov/Publications/3185.html>



90% of all out-of-hospital sudden cardiac arrest events in the US are fatal.
(American Heart Association)

What is an AED?

An AED, or automated external defibrillator, is used to help those experiencing sudden cardiac arrest. It's a sophisticated, yet easy-to-use medical device that can analyze the heart's rhythm and if necessary, deliver an electrical shock, or defibrillation, to help the heart re-establish an effective rhythm.

Why should I consider an AED?

- Cardiovascular disease is the leading global cause of death¹
- More than 356,000 people will suffer an out-of-hospital cardiac event this year¹
- The average response time for first responders once 911 is called is 8-12 minutes²
- Odds of survival drop by 10% with each minute that passes without an AED²

1. American Heart Association. 2021. CPR Facts & Stats.
https://www.heart.org/idc/groups/ahamh-public/@wcm/@sop/@smd/documents/downloadable/ucm_470704.pdf

2. American Red Cross. 2021. AED Training.
<https://www.redcross.org/take-a-class/aed/aed-training>




What is Sudden Cardiac Arrest?

Sudden cardiac arrest (SCA) often occurs without warning, triggered by an electrical malfunction that causes an irregular heartbeat (arrhythmia). With its pumping action disrupted, the heart can't pump blood to the brain, lungs and other organs. Seconds later, a person loses consciousness and has no pulse and death occurs within minutes if the victim doesn't receive treatment.¹

Is it the same as a heart attack?

People often use these terms interchangeably, but they are different. A heart attack occurs when blood flow to the heart is blocked. Sudden cardiac arrest occurs when the heart malfunctions and suddenly stops beating. A heart attack is a "circulation" problem and sudden cardiac arrest is an "electrical" problem.¹

Cardiac Arrest **vs.** Heart Attack¹



1. American Heart Association. 2021. Cardiac Arrest vs. Heart Attack. <https://cpr.heart.org/en/resources/cardiac-arrest-vs-heart-attack>

What should I do if someone at work experiences SCA?

- **Call 911** for emergency medical services
- **Locate an AED** if one is available and use it as soon as it arrives
- **Begin CPR immediately** and continue until professional medical help arrives

By performing Hands-Only CPR, you can still double or even triple a victim's chance of survival.*

*American Heart Association

AED Service Agreement

Setup in 3 easy steps:

1 PLANNING

- You determine the type, location and number of AEDs needed for your business

2 IMPLEMENTATION


- Your AEDs and accessories you select are professionally installed

3 MAINTENANCE

- Pad and battery replacement
- Software upgrades
- Ongoing LifeREADY 360™ support for medical direction, tracking, registration with EMS and post-event guidance*

Your dedicated Cintas professional will perform inspections of pad and battery expiration dates

* Customer must obtain a prescription for the AED Device(s) and the AED Device(s) must be subject to medical oversight. Cintas is not a health-care provider and does not issue the prescription, provide medical oversight, or provide other related health-care or medical services. Rather, when executing a Service Agreement, Customer receives access to a third-party, at no additional cost to Customer, who provides a prescription for the AED Device(s), provides medical oversight you need it most.



Program Management

When implementing an AED program*, businesses might consider:




Pinpointing places where people tend to be close together



Using a two-minute response time as a guideline



Identifying possible obstacles to AEDs



Making sure the AED is easily visible and accessible



Ensuring your AED is inspected on a regular basis



Updating emergency response plans and training your staff

* Newman, Christenson, Eisenberg. Community CPR-AED Programs.
Retrieved from <http://www.sca-aware.org/community-cpr-aed-programs>

* Cintas expressly disclaims any responsibility or obligation to provide any risk or hazard analysis of any kind or type, including (but not limited to) recommendations regarding the type, number, or placement of AEDs at Customer's facility. Customer is solely responsible for choosing the type and number of AEDs needed and the placement of those AEDs within its facility(ies) and for notifying Cintas of any faults, alarms, or indications that the AED is not functioning properly in between Cintas's periodic service visits.



LifeREADY 360™ powered by RescueStat

Your subscription to LifeREADY 360 includes a medical prescription along with the following:

- **Medical Direction or Physician Oversight** — by using the LifeREADY 360 portal, your AED program and activity will be overseen by the professional RescueStat team.
- **AED Placement Reporting/EMS Registrations** — your AED location information will be registered with local and/or state EMS agencies (as may be required by law) in an effort to enhance and expedite emergency responder response times to a cardiac emergency.
- **Tracking Software** — LifeREADY 360 empowers you to review your ongoing maintenance and AED readiness, by logging your AED checks, expiration dates and even tracking certified first responders.
- **AED Program Design Guide** — LifeREADY 360 provides a “design guide” you can use to develop a program that meets your needs.

The screenshot shows the LifeREADY 360 Customer Portal. The browser address bar displays `safety.rescuestat.com/CustomerPortal/EquipmentOverview`. The page header includes the LifeREADY 360 logo, a phone number (833) 507-6877, a 'REPORT EVENT' button, and a 'My Account' link with a 'Logout' option. A left sidebar contains a 'Customer' profile icon and navigation links for Home, Dashboard, Locations, Reports, Documents, and AED Laws. The main content area features a header for 'XYZ Co. - LR360_ID' with an address (6800 Cintas Blvd, Mason, Ohio 45040) and a renewal date (8/1/2023). Below this are four status cards: 'INSPECTIONS' (green checkmark), 'EQUIPMENT' (green checkmark), 'TRAINING' (yellow warning triangle), and 'REMOTE MONITORING' (N/A). A key indicates green for 'Current', yellow for 'Attention', and red for 'Expired'. At the bottom, there is an 'Equipment Overview' section with a search bar and a table with columns: AED Image, AED Make Model, Serial Number, Placement Description, Type, Expires, and Status. The table shows 10 entries.

Customer must obtain a prescription for the AED Device(s) and the AED Device(s) must be subject to medical oversight. Cintas is not a health-care provider and does not issue the prescription, provide medical oversight, or provide other related health-care or medical services. Rather, when executing a Service Agreement, Customer receives access to a third-party, at no additional cost to Customer, who provides a prescription for the AED Device(s), provides medical oversight you need it most.

Why Cintas?

- Nearly 800 lives saved with Cintas-deployed AEDs since 2007
- 80 locations across North America with over 1,600 dedicated service professionals
- The largest van-delivered first aid supplier in North America
- One of the largest providers of the American Heart Association First Aid and CPR/AED course, certifying over 300,000 individuals every year
- More than 200 certified instructors
- Over 350 safety training courses to choose from
- Exclusive relationships with industry-leading AED manufacturers



Consider Cintas a resource to help you implement an AED program — so you can be prepared when it matters most.



CINTAS RECOGNIZES BUSINESSES
WHO INVEST IN AEDS
AND HAVE THE **CONFIDENCE**
TO USE THEM TO **SAVE A LIFE**

The AED 3[®] by ZOLL[®]

TOUCH SCREEN DISPLAY

- Full-color rescue images with text
- One-of-a-kind touch screen for quick access to AED information

WIFI CONNECTED

- Cloud connectivity can enable automatic reporting of device status
- Quickly transmit data to medical professionals

ENHANCED REAL CPR HELP




- Real CPR Help technology helps support rescuer with verbal prompts
- CPR cycle timer and compression measurements
- User will be alerted if shock is needed

SMART BATTERY

- Power level always shown on screen

INTEGRATED PEDIATRIC RESCUE

- Uni-pads applicable to both children and adults
- To treat a child, simply activate Child Mode



50 states in the U.S. have Good Samaritan laws in place to protect bystanders who use an AED to resuscitate an SCA victim.

Recreation Law. May 2014. Good Samaritan Laws by State. <https://recreation-law.com/2014/05/28/good-samaritan-laws-by-state/>

Accessories



ELECTRODE PAD

For use with adults and children, simply change pad placement and press the pediatric mode button

616852 CPR Uni-padz™ electrodes



WALL STORAGE


Keeps AED visible and secure

19140 XL AED Cabinet w/ Alarm (20" x 20" x 12")

19150 Deluxe AED Cabinet w/ Alarm (17.5" x 17.5" x 9")

616854 ZOLL® AED 3® Unit Only Wall-mount

616855 ZOLL® AED 3® Carry Case Wall-mount



AED SIGNAGE

Informs employees and visitors of the availability and location of AEDs

341810 6" x 9" 3D-sign

341790 10" x 14" flat sign



REPLACEMENT BATTERY

Standby life of five years (installed in the unit)

616853 ZOLL® AED 3® Disposable Lithium Manganese Dioxide Battery



CARRY CASES

Transport your ZOLL® AED 3®, an extra set of pads, extra battery and other accessories during emergencies

616856 Standard Carry Case

(holds extra set of pads)

616857 Large Hard-shell Carry Case

(holds extra set of pads and one spare battery)

616858 Slim Hard-shell Carry Case

(holds one spare battery)

The AED Plus[®] by ZOLL[®]



RESCUE WITH CONFIDENCE

- Real CPR Help technology measure compression depth for CPR accuracy
- Provides verbal cues like "Push Harder" or "Good compressions"

ACCURATE PAD PLACEMENT

- The CPR D-pads feature a unique one-piece design to ensure accurate pad placement

LONG LASTING AND DURABLE

- Designed to withstand harsh environments and extreme temperatures
- Highly resistant to dust and moisture
- Industry-leading batteries and pads that last five years



When CPR and an AED are administered to a victim within the first two minutes of a cardiac event, the chances of survival increase 90%.

Safety and Health Magazine. May 2016. Automated external defibrillators at work. <https://www.safetyandhealthmagazine.com/articles/14041-automated-external-defibrillators-at-work>

Accessories



ELECTRODE PADS

For adults and children 8 years or older or over 55 lbs.
8306150 CPR-D-padz® electrodes For infants and children 8 years or under or under 55 lbs.
8306188 Pedi-padz® II pediatric electrodes




WALL STORAGE

Keeps AED visible and secure
19140 XL AED Cabinet w/ Alarm (20" x 20" x 12")
19150 Deluxe AED Cabinet w/ Alarm (17.5" x 17.5" x 9")
577129 ZOLL® AED Plus® Wall-mount Bracket



CARRY CASES

Transport your ZOLL® AED Plus®, an extra set of pads, extra battery and other accessories during emergencies
ZM8000080201 Black Carry Bag (holds extra set of electrodes and a sleeve of batteries)



AED SIGNAGE

Informs employees and visitors of the availability and location of AEDs
341810 6" x 9" 3D-sign
341790 10" x 14" flat sign



REPLACEMENT BATTERY

Standby life of up to five years (installed in the unit)
8306151 Duracell 123A Lithium Batteries (Sleeve of 10)

The Lifeline VIEW by Defibtech



INDICATOR LIGHT

- Daily self-test to ensure functionality
- Blinking red indicates low battery
- Flashing green indicates ready status

INTERACTIVE VIDEO DISPLAY

- Step-by-step videos and text prompts
- Interactive display coaches user through defibrillation
- Provides guidance on how to apply pads and perform CPR

VOICE PROMPTS

- Verbal cues that lead user through the process of a rescue

STATUS SCREEN

- One-Touch allows you to check unit status without powering on
- Instant verification of expiration dates

defibtech

AED status	OK
Battery status	Expires 02/2025
Pads status	Adult Expires 01/2022
AED S/N	4300xxxxx
Battery S/N	4300xxxxx
Software version	Vx.x



Waiting for the arrival of emergency medical system personnel results in only 5-7% survival.

OSHA. May 2003. Saving Sudden Cardiac Arrest Victims in the Workplace. <https://www.osha.gov/sites/default/files/publications/osh3185.pdf>

Accessories



ELECTRODE PADS

For adults and children 8 years or older or over 55 lbs.

DEDDP2001 Lifeline VIEW Adult Pads

For infants and children under 8 years old

DEDDP2002 Lifeline VIEW Pediatric Pads



WALL STORAGE

Keeps AED visible and secure

19110 Cintas Basic AED Cabinet w/Alarm (18" x 17.25" x 4")

19120 Lifeline AED Cabinet (No Alarm) (18" x 17.25" x 4")

132020 Wall-mount bracket (for carrying cases)




EMERGENCY RESPONSE SUPPLIES

Easily transport first aid supplies with your AED in the event of an emergency

21859 Lifeline Prep & Response Kit (absorbent towel, protective gloves, CPR mouth barrier, antiseptic wipe, medical scissors and disposable razor in carry case)

21880 Lifeline Prep & Response Kit Refill

21790 Grab-N-Run Quick Response Kit w/ Trauma Fill



AED SIGNAGE

Informs employees and visitors of the availability and location of AEDs

341810 6" x 9" 3D-sign

341790 10" x 14" flat sign



REPLACEMENT BATTERY

Standby life of four years (installed in the unit)

DEDCF2003 4-Year Lifetime

Lifeline VIEW Battery Pack



CARRYING CASE

Holds Lifeline VIEW, extra set of pads, extra battery and other accessories

21780 Ballistic Nylon Soft Carrying Case

The Lifeline by Defibtech



TWO-BUTTON OPERATION

TEXT, LIGHTS AND VOICE GUIDES

- Voice-guided prompts help the responder step-by-step
- Visual display of progress using lights

BUILT-IN METRONOME

- Produces a steady beat during CPR

INDICATOR LIGHT

- Daily self-test to ensure functionality
- Blinking red indicates low battery
- Flashing green indicates ready status

DATA CARD SUPPORT

- Data Card technology makes for easy upgrades on-site
- Downloads data for event review when deployed

DEFIBTECH Lifeline AED DEFIBRILLATOR

check pads
do not touch patient
analyzing

SHOCK

defibtech DDC-2001 Data Card

SHOCK

SHOCK



About 400 workplace deaths from cardiac arrest are reported to OSHA each year.

OSHA, May 2003. Saving Sudden Cardiac Arrest Victims in the Workplace. <https://www.osha.gov/sites/default/files/publications/osh3185.pdf>

Accessories



ELECTRODE PADS

For adults and children 8 years or older or over 55 lbs.
DEDDP100P1 Lifeline Adult Pads

For infants and children under 8 years old
DEDDP200P Lifeline Pediatric Pads



WALL STORAGE

Keeps AED visible and secure

19110 Cintas Basic AED Cabinet w/Alarm (18" x 17.25" x 4")

19120 Lifeline AED Cabinet (No Alarm) (18" x 17.25" x 4")

132020 Wall-mount bracket (for carrying cases)




EMERGENCY RESPONSE SUPPLIES

Easily transport first aid supplies with your AED in the event of an emergency

21859 Lifeline Prep & Response Kit (absorbent towel, protective gloves, CPR mouth barrier, antiseptic wipe, medical scissors and disposable razor in carry case)

21880 Lifeline Prep & Response Kit Refill

21810 Grab-N-Run Quick Response Kit w/ Trauma Fill



AED SIGNAGE

Informs employees and visitors of the availability and location of AEDs

341810 6" x 9" 3D-sign

341790 10" x 14" flat sign



REPLACEMENT BATTERY

Standby life of four years (installed in the unit)

DEDCF200 5-Year Lifetime Lifeline Battery Pack



CARRYING CASE

Holds Lifeline, extra set of pads, extra battery and other accessories

DEDAC100 Ballistic Nylon Soft Carrying Case

American Heart Association® Heartsaver® First Aid, CPR/AED Training



This course includes hands-on instruction with engaging training your employees will enjoy and physical documentation for all attendees. Schedule at your convenience including night and weekend classes.

Also available as an online course with in-person CPR skills check.

LEARN HOW TO:

- Recognize and treat adult emergencies prior to EMS arrival
- Administer basic First Aid
- Become certified to perform adult CPR
- Effectively use an AED to deliver a life-saving shock
- Pediatric CPR module available

Federal Regulation:

29 CFR 1910.151(b) "In the absence of an infirmary, clinic, or hospital in the near proximity to the workplace, which is used for the treatment of all injured employees, a person or persons shall be adequately trained to render First Aid."

Course #650011

Depending on the nature and location of the worksite, federal, state, and local regulations may require that employers offer specific training to their employees. Employers or businesses have sole responsibility to examine their operations to determine the type of training that is required or appropriate. Cintas does not provide recommendations or guidance regarding which safety and training programs are required or appropriate. The services are provided for educational purposes only.

Emergency Oxygen

Providing supplemental oxygen early in the course of a heart attack or stroke can increase the amount of oxygen reaching the heart and brain, and in turn help increase chances of survival and quick recovery. The average response time for EMS to arrive on the scene is 8-12 minutes.²

The R15™ by Rapid Oxygen

- Provides 15 minutes of humidified oxygen at 6L/min
- Easy to use, no maintenance required
- No prescription required
- No electricity or batteries required
- Non-explosive and requires no specialized storage
- Safe to store and use in all public areas



617075 - The R15™
617077 - The R15™ Shelving Unit

Easy to use – activates in just 3 simple steps



1. https://professional.heart.org/idc/groups/ahamam-public/@wcm/@sop/@smd/documents/downloadable/ucm_505473.pdf
2. <https://medicalnewsbulletin.com/response-time-emergency-medical-services/>

When it comes to the heart, there's no room for doubt.
We've featured some of the industry's most innovative AEDs
here, but additional accessories and units may be available.

To learn more, visit **cintas.com** or call 877.973.2811

cintas.com | 877.973.2811

UNIFORMS | FACILITY SERVICES | FIRST AID & SAFETY | FIRE PROTECTION

FAS-AED-PS rev.4/22

